

## Admitting

1. Use a barcode reader to scan the patient's barcode, or Touch **MENU** → **PATIENT LIST** → **NEW PATIENT**.
2. Then touch each desired data field to enter appropriate patient information.
  - Remember to select **ENT** or **SET** after each selection to save the selected data.




## Saving and Sending to EMR

1. After collecting the desired vital signs of the admitted patient, Touch the **SAVE** key → select/unselect parameters to send (unselected items will turn grey) → **SAVE & SEND** key.
  - Remember, in spot mode, the patient will be discharged automatically after you SAVE the information. In continuous mode, you must discharge (QUIT) the patient before moving to the next patient.



## Discharging

1. Touch **MENU** → **PATIENT INFO** → **QUIT PATIENT**.
2. Touch the **YES** key on the confirmation window.
  - The data is then cleared on the Home screen.



## Reviewing Vital Signs (V/S)


### Spot Table Review

1. Touch **MENU** → **SPOT TABLE**.

### Trend Table Review

2. Touch **MENU** → **TREND TABLE**.

Review Windows



## Spot Table Window

When NIBP or TEMP measurement is completed or when (Save) is pressed, the measuring values may be displayed on the Spot Table screen.

1. Patient information filter area  
Tap the dialog box to choose the patient name or ID. Tap **ALL** to display all data.


2. Send key  
Tap **SEND** to send the record you choose.

3. Delete key  
Tap **DELETE** to delete the record you choose.

4. Detail icon  
Tap to display the Spot Data Detail window.

5. Record area  
The records which have been sent are in grey. The records which have not been sent are in black.

6. Check box  
Tap the check box to choose the record you need.



## Trend Table Window (Continuous Mode Only)

1. Patient information filter area  
Tap the dialog box to choose the patient name or ID. Tap **ALL** to display all data.

2. Time interval filter area  
Tap the dialog box to choose 1 min, 5 min, 10 min, 15 min, 30 min or 1 h.

3. Record area  
The records which have been sent are in grey. The records which have not been sent are in black.

