



# ALARM MANAGEMENT & BEST PRACTICES

# History of Alarms



The word alarm originates from the Latin 'ad arma' or French 'a l'arme' which translates to 'to your weapons'



Some of the first documented alarms dates back to the Middle Ages



An alarm is a warning that results from a measurement variance and indicates a deviation from normal state



Alarms are essential and necessary in the monitoring of vital signs necessary to support life



The word 'alarm' indicates a call for immediate action or defense



Alarms are intended to prevent patient harm by providing rapid reaction to critical situations **BUT** only if they are not false alarms

# ALARM FATIGUE

Clinicians become desensitized, overwhelmed or immune to the sound of an alarm

## What is Alarm Fatigue?



Alarm fatigue is the direct result of the constant bells, blips and alarm signals emitted by medical devices



Care givers may become 'immune' to these sounds which increases the risk of these alarms being absorbed into the auditory landscape of hospital corridors and subsequently being ignored<sup>2</sup>

## Fatigued Clinicians May:



**Turn down  
alarm volume**



**Turn off  
alarm**



**Adjust alarm  
settings**

**These actions can have serious or fatal consequences**

# Alarm Fatigue



## Hazard Faced by Hospitals

“In a hospital setting, the frequency of alarms poses a risk of some hospital staff becoming **desensitized** to the constant beeps and in the worst cases, **lowering the volume too much**”

- ECRI Institute

# Patient Safety Issue



**Sentinel Alert:** (OAKBROOK TERRACE, Ill. – April 8, 2013)  
The constant beeping of alarms and an overabundance of information transmitted by medical devices such as ventilators, blood pressure monitors and ECG (electrocardiogram) machines is creating “alarm fatigue” that puts hospital patients at serious risk, according to a [Sentinel Event Alert](#) issued Today by The Joint Commission.



**#2**  
**Missed Alarms**  
**Technology Hazard**

## Missed Alarms Can Have Fatal Consequences!

# The Joint Connection

The Scope of the Problem

## Medical Device Alarm Safety

**100s**

of alarm signals per patient per day

**1,000s**

of alarm signals on each unit

**10,000s**

of alarm signals throughout a hospital per day



**85-99%**

**of alarm signals** don't require clinical intervention

# The Joint Connection

## Patient Events



# National Patient Safety Goals on Alarms

The Joint Commission Announces 2014 National Patient Safety Goal In June 2013

## PHASE 1

**(2014 to January 1, 2016)**

---

Establish alarm system safety as hospital priority

Identify most important alarm system to manage (EC.02.04.01)

## PHASE 2

**(began January 1, 2016)**

---

Establish policies and procedures for managing the alarms identified

Educate staff and communicate changes



# NPSG on Alarm Management – Phase I

NPSG on Alarm Management In Phase I (beginning January 2014)



Establish alarm management as an organizational priority



Identify most important alarms to manage based on internal situations



Seek input from medical staff and clinical departments



Identify risks to patients due to lack of response and/or malfunction



Identify actionable alarms vs alarms contributing to noise/fatigue



Publish best practices/guidelines



**Pro-actively think:**  
Identify potential for patient harm based on internal incident history

# NPSG on Alarm Management: Phase 2

NPSG on Alarm Management In Phase II (beginning January 2016)

**Hospitals were expected to develop and implement specific components of policies and procedures that address at minimum:**



Clinically appropriate settings



Process for monitoring alarms and expectations moving forward



When alarms can be disabled



Checking individual alarm signals for accurate settings, proper operation and detectability



When parameters can be changed



Educate those in the organization about alarm policies



Who can set, who can change parameters, and who can set to "off"



Clinical Leadership responsible ensuring accuracy, safety, education and communication around new alarm policies

# Common Cause for Nuisance Alarms



Improper lead placement



Many parameters set to always sound an alarm



No schedule to change leads



Multiple **unverified alarms** being routed to nurse pager or phone



Frayed or malfunctioning lead wires



Lack of patient education, thus decreased patient compliance with wearing the device



Alarms not customized



Alarm defaults still per manufacturer recommendations, not customized to reflect unit/hospital policies



# NPSG on Alarm Management: Taking Action



**Create** or participate in an Alarm Management Committee



**Change** transmitter batteries daily or when the system indicates the need for them



**Stress** loops on ECG lead wires for tremulous patients



**Customize** individual patient alarms based on assessment and condition to assure alarms are valid



**Establish** organizational alarm defaults by unit





**Customize** alarms based on patient assessment and report/verify at shift changes



**Change** electrodes daily using recommended procedures. Add this to daily routine tasks such as baths, PM care to assure practice change

# Stress Loops for ECG Leads



# Best Practices: Alarm Settings

## Alarm settings, limits and delays



Establish appropriate default settings for hospital unit and patient population



Small changes can yield big results:  
Example: Decreasing SpO2 lower threshold by one point from 90% - 89%



Turn off **duplicate alarms**



Consider using **alarm signal delays** to all for alarm autocorrection



Ensure alarm priority is set to **actionable levels**



Consider using secondary alarm notification to improve audibility



Review high/low settings and other limits



Consider alarm escalation to increase priority

# Best Practices: Alarm Settings



Establish alarm limits and defaults based on **population served**  
(Adult ICU vs NICU vs Emergency Department)



Create a process to customize alarm settings based on **individual patients**



One size **does not fill all**

## Staff Education



Educate clinicians on their roles related to alarm management



Empower staff to manage nonactionable alarms by changing limits to actionable levels (in accordance with organizational policy)



Ensure staff are trained and competent in recognizing and troubleshooting equipment alarm signals



Encourage staff to review trend data for repetitive alarms, especially during sleep. A perceived false alarm may be a sleep apnea patient with multiple clinical relevant alarms that self correct when the nurse enters the room and wakes the patient

# Best Practices: Alarm Settings



## Educate

the patients and families about the physiologic monitoring systems and their role in patients safety and alarm malmanagement



## Encourage

patients and families to notify staff when an alarm signal is not being addressed in a timely manner




## Consult

with other hospitals to determine where they have set their default physiologic monitor settings



# Best Practices



## Waveform Artifact/ECG/SPO2



Review proper skin and electrode prep



Maintain regular schedule for changing electrodes



Inspect re-usable lead wires to assure they are intact, proper connections and frayed wires



Consider use of disposable lead wires and SPO2 probes



Check sensor placemat for adhesion



Stress looping

# Smoothing & Time Delays

A large percentage of clinical alarms are caused by only a mild threshold violation



Implementation of time delays on certain parameters have proved to be effective in alarm reduction



Research has shown that a 15 second smoothing algorithm on SpO<sub>2</sub> reduced false alarms by 50%



Smoothing Algorithms - Remove artifact and smooth the vital sign data collected



Program average periods for heart rate/pulse rate (usually 0-10 seconds), SpO<sub>2</sub> and respiration rates (usually 0-30 seconds)

# Responding to Notification Systems



**Intervene to the Notification in a Timely Manner**



**Assess the PATIENT When Alarms Occur**



**Take Action to Correct Patient Problem**



**Replace probe(s)**



**Replace/Change Electrodes**



**Adjust Alarm Parameters**

Based on Clinical Assessment and Stability



**Assess Trends Proactively**




**Educate Patient**

On need for continuous monitoring to increase patient understanding and nursing's commitment to patient safety

# How Do I Measure My Changes?

Data – What to Look for in Alarm Reports



# How Do I Measure My Changes?

Data – What to Look for in Alarm Reports



## Real-time Data Or Metrics For Dashboard Reporting



Typically Most Useful to Point-of Care Nurses



Shifts with Most Alarms



Alarms by Nursing by Unit



Time of Week for Most Alarms



Alarms by Bed, Shift or Time of Day



Type of Alarm Sounding Most Often

# How Do I Measure My Changes?

How Does Committee Determine Actionable vs. Non Actionable Alarms?

## Critical Factors



Include point of care professionals



Consider how information about alarms and types of alarms could affect work environments



<b>Staffing</b>	<b>Workflow Analysis</b>	<b>Alarms Settings and Management</b>
<b>Immediate Impact of Alarm Flood</b> (10 or more alarms in 10 minutes)	<b>Protocol Management and Review</b>	<b>Metric Comparison Against Peers</b>

# Monitoring for Outcomes

Clinically Relevant Process to Ensure No Negative Outcomes



**Monitor Noise Level  
in the Unit**



**Survey Patient  
Satisfaction**



**Evaluate Incident Reports**

Determine if issues are a result of  
alarms management changes



**Convene Regular  
Safety Huddle**



**Monitor for Any Increase**

Increase in rapid response calls, codes,  
unplanned transfers to critical care



**THANK YOU FOR YOUR TIME!**

To download this presentation, visit [nkuniversity.org/resources](https://nkuniversity.org/resources)



# References

- <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC3672485/>
- <http://medcitynews.com/2012/03/alarm-fatigue-becomes-focus-of-fda-attention/>
- <http://medcitynews.com/2012/03/alarm-fatigue-becomes-focus-of-fda-attention/>
- [http://www.jointcommission.org/sea\\_issue\\_50/](http://www.jointcommission.org/sea_issue_50/).
- [http://www.jointcommission.org/sea\\_issue\\_50/](http://www.jointcommission.org/sea_issue_50/).
- [http://www.jointcommission.org/sea\\_issue\\_50/](http://www.jointcommission.org/sea_issue_50/).
- Clochesy JM, Cifani L, Howe K et al. Electrode site preparation techniques: a follow-up study. Heart Lung. 1991;20:27-30
- <http://www.aacn.org/wd/practice/docs/practicealerts/alarm-management-practice-alert.pdf>
- Cospers, P., Zellinger, A., Jacques, L., Razzano, L., Flack, M., Improving Clinical Alarm Management: Guidance and Strategies. BI&T 2016;51,109-115
- Cospers, P., Zellinger, A., Jacques, L., Razzano, L., Flack, M., Improving Clinical Alarm Management: Guidance and Strategies. BI&T 2016;51,109-115
- Cospers, P., Zellinger, A., Jacques, L., Razzano, L., Flack, M., Improving Clinical Alarm Management: Guidance and Strategies. BI&T 2016;51,109-115
- <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC3672485//>
- Cospers, P., Zellinger, A., Jacques, L., Razzano, L., Flack, M., Improving Clinical Alarm Management: Guidance and Strategies. BI&T 2016;51,109-115
- Cospers, P., Zellinger, A., Jacques, L., Razzano, L., Flack, M., Improving Clinical Alarm Management: Guidance and Strategies. BI&T 2016;51,109-115
- Cospers, P., Zellinger, A., Jacques, L., Razzano, L., Flack, M., Improving Clinical Alarm Management: Guidance and Strategies. BI&T 2016;51,109-115
- Cospers, P., Zellinger, A., Jacques, L., Razzano, L., Flack, M., Improving Clinical Alarm Management: Guidance and Strategies. BI&T 2016;51,109-115